



  **Buchanan Orthotics**

# Steprite Footwear

LADIES

GENTS

KIDS



# Free Your Feet

Buchanan Orthotics is the largest Orthotic contractor in Scotland, supplying both the NHS and private patients with a full range of orthotic products. A long established bespoke footwear and orthoses manufacturer, utilised this expertise to develop our Steprite range of modular footwear which is designed to meet the demand for quality footwear.

Since its introduction Steprite footwear has been very successful, and our designers are constantly updating colours and introducing new styles to allow more patient choice. A greater choice of new styles, colours and leathers now makes owning a pair Steprite shoes more of a fashion statement, whilst still accomodating patients requirements and expectations.





## CONTENTS



Page 3 Ladies Shoes

Page 8 Ladies Low Opening Shoes

Page 10 Ladies Stretchable Shoes

Page 12 Ladies Boots

Page 13 Ladies Diabetic Specifications

Page 14 Ladies Trainers

Page 17 Gents Shoes

Page 20 Gents Low Opening Shoes

Page 23 Gents Stretchable Shoes

Page 24 Gents Diabetic Specifications

Page 25 Gents Boots

Page 27 Gents Trainer

Page 28 Stock Footwear s

Page 30 Children's KP Athletic & School Shoe

Page 32 Bespoke Footwear Style Options For Ladies & Gents

Page 34 TCIs, Plastic Orthoses & FFOs

Page 36 Fabric Supports

Page 37 Colours

Page 38 Ladies Last Accommodation Charts

Page 40 Gents Last Accommodation Charts

Page 42 Ladies Contour Lasts

Page 43 Gents Contour Lasts

Page 42 Last Adaptions

Page 43 Sole Options

Page 45 Conditions of Sale





## STEPRITE FOR LADIES

Ladies Steprite footwear offers a wide range of both traditional and casual styles. As an established bespoke footwear and orthoses manufacturer we have incorporated an extensive range of fittings and colour choices to this new range, all of which are made to the highest standard of craftsmanship.

### Ailsa

Five eyelet casual with extended punched vamp panels, tongue attached externally to avoid internal seams, bambus quarter, vamp and tongue lining with padded collar and suede heel grip, lightweight through sole.

Colours	All leathers and numbuck
Lasts	Regular, Donna, Diana, Denise and Contour
Sole Colour	Black or brown
Fastening	Lace or Touch & Close Return



LACE



TOUCH & CLOSE RETURN



### Iona

Ladies casual soft leather uppers with padded tongue, Extended vamp with cord stitch design rows. External speed lace loops. Bambus vamp lining, leather quarter lining, padded collars. Lightweight through sole

Colours	All colours
Lasts	Regular, Donna, Diana, Denise and Contour
Sole Colour	Black or brown
Fastening	Lace or Touch & Close Return



LACE



TOUCH & CLOSE RETURN

## Maryland

Modern casual: Extended vamp with multi-stitch design rows. Bambus vamp lining, leather quarter lining with padded achilles and tongue, lightweight through sole.

Colours	All colours
Lasts	Regular, Donna, Diana, Denise and Contour
Sole Colour	Black or Brown
Fastening	Lace or Touch & Close Return



LACE



TOUCH &amp; CLOSE RETURN

## Sonoma

Stitch design with apron mudguard, leather quarter linings, bambus vamp linings, padded collar.

Colours	All Leathers
Lasts	Regular, Donna, Diana and Denise
Sole Colour	Black, Brown and Taupe
Fastening	Lace or Touch & Close Return



LACE



TOUCH &amp; CLOSE RETURN

### Skye

Ladies Gibson with soft leather and tweed fabric apron. Padded collars with leather quarter linings, bambus vamp linings. Lightweight sole/heel.

Colours Black, Brown, Taupe

Lasts Regular, Donna, Diana and Denise

Fastening Lace or Touch & Close



LACE



TOUCH & CLOSE

### Louise

Single strap shoe with padded collar. Lightweight sole and heel or trainer sole.

Colours All Colours

Lasts Regular, Donna and Denise

Fastening Touch & Close only



TOUCH & CLOSE



TOUCH & CLOSE

### Manhattan

Multi colour patchwork apron with mudguard, leather quarter linings, Bambus vamp linings.

Colours	All Leathers
Lasts	Regular, Donna, Diana, Denise and Contour
Sole Colour	Black, Brown, Taupe
Fastening	Lace or Touch & Close Return



LACE



TOUCH & CLOSE RETURN

### Phoenix

Plain vamp and quarter, bambus vamp lining, leather quarter lining with padded collars.

Colours	All Leathers, Suede, Nubuck
Lasts	Regular, Donna, Diana, Denise and Contour
Sole Colour	Black, Brown, Taupe
Fastening	Lace or Touch & Close Return



LACE



TOUCH & CLOSE RETURN



## Orlando

Seam free fronted shoe with punch design, bambus vamp lining, leather quarter linings with padded collar.

Colours	All Leathers and Suedes
Lasts	Regular, Donna, Diana, Denise and Contour
Sole Colour	Black, Brown, Taupe
Fastening	Lace or Touch & Close Return



LACE



TOUCH &amp; CLOSE RETURN

## Illinois shoe & sandal

Casual style low opening shoe with touch and close. Lightweight through sole.

Colours	All Leathers and Suedes
Lasts	Regular Donna, Diana and Denise
Sole Colour	Black or Brown
Fastening	Touch & Close



SHOE



SANDAL

## Low Opening Styles

## Boston

Low opening shoe. Casual style shoe with apron design rows and contrasting side panel. Leather quarter linings, bambus vamp linings and padded collar.

Colours	All Leathers, Suedes and Nubucks
Lasts	Regular, Donna, Diana, Denise and Contour
Sole Colour	Black, Brown, Taupe
Fastening	Lace or Touch & Close Return



LACE



TOUCH &amp; CLOSE RETURN

## Alabama

Very low opening to cap for easy entry, soft leather uppers, bambus tongue and vamp linings, sueded counter lining, seam free vamp lining.

Colours	All Leathers
Lasts	Regular, Donna, Diana, Denise and Contour
Sole Colour	Black, Brown, Taupe
Fastening	Lace or Touch & Close Return



LACE



TOUCH &amp; CLOSE RETURN

## Arran

Ladies low opening casual shoe with cord design row. Lace retaining loops to prevent tongue migration. Bambus seamless vamp lining with sueded heel grip and padded collars. Softee leather.

Colours	All Leathers
Lasts	Regular, Donna, Diana, Denise and Contour
Sole Colour	Black, Brown, Taupe
Fastening	Lace or Touch & Close Return



LACE



TOUCH &amp; CLOSE RETURN

## Phoenix Low

Low opening shoe, approx 20mm lower than standard. Plain vamp and quarter leather quarter lining, bambus vamp lining with padded collar.

Colours	All Leathers and Nubucks
Lasts	Regular, Donna, Diana, Denise and Contour
Sole Colour	Black, Brown, Taupe
Fastening	Lace or Touch & Close Return



LACE



TOUCH &amp; CLOSE RETURN

# Stretchable Styles

Stretch leather only available in Black, Brown , Navy, Taupe and Leopard Skin

## Jura

Ladies side panel stretch leather shoe. Linings: polycotton, suede heel grip, padded collar. Lightweight through sole with skived and wedged heel for better shock absorption.

Colours	Black, Brown, Navy and Taupe
Lasts	Regular, Donna, Diana ,Denise and contour
Fastening	Lace or Touch & Close Return



LACE



TOUCH & CLOSE RETURN



## Portree

Ladies gibbon with stretch leather or neoprene apron, leather quarter, counter linings, polycotton, vamp lining with padded collar.

Colours	Black, Brown, Navy and Taupe
Lasts	Regular, Donna, Diana , Denise and Contour
Sole Colour	Black or Brown
Fastening	Lace or Touch & Close Return



LACE



TOUCH & CLOSE RETURN

### Vermont

Ladies style with plain neoprene vamp. Polycotton with foam cushioning vamp lining, leather quarter lining and padded collars. Available with or without toe puff.

Colours	Neoprene: Black, Brown, Navy, Burgundy, Light Beige
Lasts	Regular, Donna, Diana, Denise and Contour
Sole Colour	Black or Brown
Fastening	Lace or Touch & Close Return



TOUCH & CLOSE



TOUCH & CLOSE RETURN



TOUCH & CLOSE

### Louise

Single strap shoe with padded collar, hard wearing through flexible sole, leather stretch side panels

Colours	Black, Brown, Navy and Taupe
Lasts	Regular, Donna, Diana and Denise
Fastening	Touch & Close

**Delaware**

Turn collar ankle boot. Extended vamp with design, stitching on leg, soft leather, bambus vamp/leg lining, leather collar lining (also available with fur collar though this is non returnable).

Colours	All Leathers
Lasts	Regular, Donna, Diana, Denise and Contour
Fastening	Lace or Touch & Close Return



LACE



TOUCH &amp; CLOSE RETURN

**Ohio**

Plain vamp with design stitch on facing bambus vamp lining, leather leg lining with padded collar. Leg height approx 130mm.

Colours	All Leathers
Lasts	Regular, Donna, Diana, Denise and Contour
Fastening	Lace or Touch & Close Return



LACE



TOUCH &amp; CLOSE RETURN

### Ohio Low

Plain vamp with design stitch on facing, bambus vamp lining, leather leg lining with padded collar. Leg height approx 130mm with approximately 20mm lower opening.

Colours	All Leathers
Lasts	Regular, Donna, Diana, Denise and Contour
Fastening	Lace or Touch & Close Return



LACE



TOUCH & CLOSE RETURN

## Diabetic Specifications

All footwear can be made to diabetic footwear specification in contour last or all regular lasts. Diabetic specification, Softee leather, Go Guard lining, Seamless Vamp Lining, Rim toe puff and padded collars.



### Panama Diabetic Concept

Casual upper design to standard diabetic specification, extended seam free vamp and through rocker sole as standard.

Colours	All Leathers
Lasts	Contour or Regular
Sole Colour	Black only
Fastening	Lace or Touch & Close Return



LACE



CONTOUR TOE SHAPE



REGULAR TOE SHAPE

**Stroller Lady**

Trainer shoe with a modern wing cap design on front and a side flash detail.

Colours	All Colours
Lasts	All Lasts
Fastening	Lace or Touch & Close Return



LACE



TOUCH &amp; CLOSE RETURN

**Austin Lady**

Modern extended wing cap trainer boot in soft leather. Leather leg lining, bambus vamp lining with padded collar. Lightweight trainer sole.

Colours	All Colours
Lasts	Ladies Regular, Donna, Diana and Denise
Sole Colour	Black or White
Fastenings	Lace or Touch & Close Return



LACE



TOUCH &amp; CLOSE RETURN



### Florida Lady

Trainer shoe with long quarter flash and stitch rows. Lightweight trainer sole. Leather quarter linings and bambus vamp and tongue lining with padded collar and tongue.

Colours	All Colours
Lasts	Ladies Regular, Donna, Diana, Denise and Contours
Fastenings	Lace or Touch & Close Return



LACE



TOUCH &amp; CLOSE RETURN

### Street Lady

Low opening casual trainer shoe.

Colours	All Colours
Lasts	All Lasts
Fastening	Lace or Touch & Close



LACE



TOUCH &amp; CLOSE



## STEPRITE FOR GENTS

The Steprite range of gents footwear is made to the highest standard and provides a wide range of styles. Each is offered in a variety of fittings and colours.



### Dakota

Stitch design apron front and facings, leather quarter linings, bambus vamp lining, padded collars.

Colours	All Leathers and Nubuck
Lasts	Regular, Albert, George, Richard and Contour
Sole Colour	Black or Brown
Fastening	Lace, Touch & Close or Touch & Close Return



TOUCH & CLOSE



TOUCH & CLOSE RETURN



### Oregon

Plain vamp shoe with punch design on facings and counters. Bambus vamp lining – leather quarter lining with padded collar.

Colours	All Leathers, Suedes, Nubuck
Lasts	Regular, Albert, George, Richard and Contour
Sole Colour	Black or Brown
Fastening	Lace, Touch & Close or Touch & Close Return



TOUCH & CLOSE



TOUCH & CLOSE RETURN



### Stuart

Plain vamp shoe with punch design on facings, counter and toe. Leather quarter lining with padded collar

Colours	Jeans
Lasts	Regular, Albert, George and Richard
Sole Colour	Black or Brown
Fastening	Lace, Touch & Close or Touch & Close Return



TOUCH & CLOSE



TOUCH & CLOSE RETURN



### Bronx

Gibson facings with cap, Oxford style vamp. Gibson opening for ease of entry and fitting. Leather quarter and tongue lining. Bambus vamp lining and padded collar.

Colours	All Leathers and Nubuck
Lasts	Regular, Albert, George, Richard and Contour
Sole Colour	Black or Brown
Fastening	Lace, Touch & Close or Touch & Close Return



TOUCH & CLOSE



TOUCH & CLOSE RETURN



## Utah

Classic plain fronted shoe with seam-free vamps and padded collars, A cushioned durable bambus vamp lining and leather quarter linings.

Colours	All Leathers and Nubuck
Lasts	Regular, Albert, George, Richard and Contour
Sole Colour	Black or Brown
Fastening	Lace, Touch & Close or Touch & Close Return



TOUCH & CLOSE



TOUCH & CLOSE RETURN



## Arizona

Classic plain vamp and quarter standard opening, leather quarter linings, bambus vamp lining, padded collars.

Colours	All Leathers and Nubuck
Lasts	Regular, Albert, George, Richard and Contour
Sole Colour	Black or Brown
Fastening	Lace, Touch & Close or Touch & Close Return



TOUCH & CLOSE



TOUCH & CLOSE RETURN

# Low Opening Styles



## Oakland Low

Low opening shoe with apron stitch design on vamp. Opens approx 20mm lower than standard. Bambus vamp lining, leather quarter lining with padded collars.

Colours	All Leathers
Lasts	Regular, Albert, George, Richard and Contour
Sole Colour	Black or Brown
Fastening	Lace, Touch & Close or Touch & Close Return



TOUCH & CLOSE THROUGH LOOP



TOUCH & CLOSE RETURN



## Jackson

Low opening, plain vamp shoe with punch and stitch detail on quarter and counter. Leather quarter lining, sueded counter lining, bambus vamp lining leather with padded collar.

Colours	All colours
Lasts	Regular, Albert, George, Richard and Contour
Sole Colour	Black or Brown
Fastening	Lace, Touch & Close or Touch & Close Return



TOUCH & CLOSE



TOUCH & CLOSE RETURN



## Campbell

Mens low opening Gibson.

Colours	All colours
Lasts	Regular, Albert, George, Richard and Contour
Sole Colour	Black or Brown
Fastening	Lace, Touch & Close or Touch & Close Return



TOUCH & CLOSE



TOUCH & CLOSE RETURN



## Michigan

Modern open to toe Gibson. Soft leather. Very low opening for easy access, seam-free vamp and vamp lining. Bambus vamp/quarter lining. Sueded counter lining. Padded collar.

Colours	All colours
Lasts	Regular, Albert, George, Richard and Contour
Sole Colour	Black or Brown
Fastening	Lace, Touch & Close or Touch & Close Return



TOUCH & CLOSE



TOUCH & CLOSE RETURN



### Arizona Low

Opens approx 20mm lower than standard opening. Classic plain vamp and quarter, leather quarter lining, bambus vamp lining, padded collar.

Colours	All Leathers and Nubuck
Lasts	Regular, Albert, George, Richard and Contour
Sole Colour	Black or Brown
Fastening	Lace, Touch & Close or Touch & Close Return



TOUCH & CLOSE



TOUCH & CLOSE RETURN



### Huron

Mens four eyelet gibson with punch design on vamp pannels, leather quarter lining, bambus vamp lining with suede heel grip and padded collar.

Colours	All leathers
Lasts	Regular, Albert, George, Richard and Contour
Sole Colour	Black or Brown
Fastening	Lace, Touch & Close or Touch & Close Return



TOUCH & CLOSE



TOUCH & CLOSE RETURN





## Tampa

Low opening Gibson with plain neoprene vamp, leather quarter linings, polycotton vamp linings and padded collars. Opens approx 20mm lower than standard opening. Now available with or without toe puff.

Colours	All colours
Lasts	All Regular and Contour
Sole Colour	Black or Brown
Fastening	Lace, Touch & Close or Touch & Close Return



TOUCH & CLOSE



TOUCH & CLOSE RETURN



## Tain

Mens Gibson with stretch leather apron or neoprene panel, leather quarter and counter lining with padded collar, polycotton vamp lining.

Colours	Black, Brown or Navy
Lasts	Regular, Albert, George, Richard and Contour
Sole Colour	Black or Brown
Fastening	Lace, Touch & Close or Touch & Close Return



TOUCH & CLOSE



TOUCH & CLOSE RETURN

All footwear can be made to diabetic footwear specification in contour last or all regular lasts. Diabetic specification, Softee leather, Go Guard lining, Seamless Vamp Lining, Rim toe puff and padded collars.



### Kentucky

Plain extended vamp design. with diabetic specification and through rocker sole as standard.

Colours	All leathers
Lasts	All mens regular lasts
Sole Colour	Black only
Fastening	Lace or Touch & Close Return



CONTOUR TOE SHAPE



TOUCH & CLOSE RETURN



### Ontario

Gibson shoe with diabetic specification and through rocker sole as standard

Colours	Black, Brown
Lasts	All mens regular lasts
Sole Colour	Black or Brown
Material	Softee Leather
Fastening	Lace or Touch & Close Return



TOUCH & CLOSE RETURN



REGULAR TOE SHAPE

**Indiana 16 Low**

Ankle boot. Plain vamp with stitch design on facing. Bambus vamp lining, leather leg lining with padded collar. Leg height approx 130mm.

Colours	All Leathers
Lasts	Regular, Albert, George, Richard and Contour
Fastening	Lace, Touch & Close or Touch & Close Return



LACE



TOUCH &amp; CLOSE RETURN

Boot legs can be made to measure on a non-return basis

**Nebraska**

Ankle boot. Plain vamp, cushioned seam-free vamp lining, leather quarter linings. Bambus vamp linings. Padded collars.

Colours	All Leathers and Nubuck
Lasts	Regular, Albert, George, Richard and Contour
Fastening	Lace, Touch & Close or Touch & Close Return



LACE



TOUCH &amp; CLOSE RETURN

### Houston

Boot with turned over apron strip and contrasting leg panels. Also available with plain vamp leather leg linings, Bambus vamp lining and deep padded collar.

Colours	All Leathers, Nubuck Brown, Taupe
Lasts	Regular, Albert, George, Richard and Contour
Sole Colour	Black or brown
Fastening	Lace, Touch & Close or Touch & Close Return



LACE



TOUCH & CLOSE RETURN

Street Men

Low opening casual trainer shoe

Colours All Colours

Lasts All Lasts

Fastening Lace or Touch & Close



LACE



TOUCH & CLOSE

Austin Men

Modern extended wing cap trainer boot in soft leather. Leather leg lining, cambrelle vamp lining with padded collar. Lightweight trainer sole.

Colours All Colours

Lasts All Lasts

Fastening Lace or Touch & Close Return



LACE



TOUCH & CLOSE RETURN

### Florida Men

Trainer shoe with long quarter flash and stitch rows. Lightweight trainer sole. Leather quarter linings and bambus vamp and tongue lining with padded collar and tongue.

Base Colour Black, white, navy, burgundy, grey, red (leather)

Flash Colour Black, white, navy, burgundy, grey, red (leather)

Sole Colour Black or white/black

Lasts Gents Regular, Albert, George and Richard

Fastening Lace or Touch & Close Return



LACE



TOUCH & CLOSE RETURN

### Stroller Men

Trainer shoe with a modern wing cap design on front and a side flash detail.

Colours All Colours

Lasts All Lasts

Fastening Lace or Touch & Close Return



LACE



TOUCH & CLOSE RETURN

## Ladies Stock Footwear



LACE



TOUCH &amp; CLOSE RETURN

## Phoenix

Colour	Black or Brown
Sizes	3 to 8 (full sizes only)
Lasts	Standard, Wide, Deep and Slim
Fastening	Lace or Touch & Close Return



LACE



TOUCH &amp; CLOSE RETURN

## Phoenix Low

Colour	Black or Brown
Sizes	3 to 8 (full sizes only)
Lasts	Standard, Wide, Deep and Slim
Fastening	Lace or Touch & Close Return

## Gents Stock Footwear



LACE



TOUCH &amp; CLOSE RETURN

## Arizona

Colours	Black or brown
Sizes	6 to 11 (full sizes only)
Lasts	Standard, Wide, Deep and Slim
Fastening	Lace or Touch & Close Return



LACE



TOUCH &amp; CLOSE RETURN

## Arizona Low

Colours	Black or Brown (shown in dark brown)
Sizes	6 to 11 (full sizes only)
Lasts	Standard, Wide, Deep and Slim
Fastening	Lace or Touch & Close Return

Express delivery subject to stock shelf availability.

# Children's KP Athletic Trainer and School Shoe

## Designed for Children's Adaptive and AFO Footwear Needs.

The KP School Shoe and Trainer, distributed by Buchanan Orthotics is a range of children's adaptive footwear. A fashionable and durable shoe that fits easily and comfortably over a brace/afo without being oversized.

For children who require orthotics to support their legs and feet.



### Distinctive Sole

Wider tread design increases direct surface contact and promotes greater medial/lateral stability.



## The KP Athletic Training Shoe

Fashionable Footwear with Innovative Design for Children's Adaptive and AFO Footwear Needs.

	Child	Adolescent
<b>Unisex UK Sizes: Grey -</b>	<b>3-13</b>	<b>1-7</b>
<b>Black -</b>	<b>10-13</b>	<b>1-8</b>

**Please advise Child or Adolescent when ordering.**

### Additional Features and Benefits

- Geometrically-graded last ensures a precise fit with accurately sized shoe.
- Unique mid-sole is a compression-molded EVA shell with a visible cut-line designed to facilitate the modification and build-up process.
- Fray-resistant lining is longer lasting and withstands the wear of AFO use.



# The KP School Shoe



**Distinctive Sole**  
 Wider tread design increases direct surface contact and promotes greater medial/lateral stability.



Fashionable Footwear with Innovative Design for Children's Adaptive and AFO Footwear Needs.

**UK Sizes: 9-13 Child, 1-10 Adolescent**



**Detachable Cut Line**  
 For ease of adaption.  
 Available on all styles of footwear.

**Additional Features and Benefits**

- Geometrically-graded last ensures a precise fit with accurately sized shoe.
- Unique mid-sole is a compression-molded EVA shell with a visible cut-line designed to facilitate the modification and build-up process.
- Fray-resistant lining is longer lasting and withstands the wear of AFO use.

## KP Accommodation Chart

Size UK	Size US	Foot Length	Width at Joint	Joint Circum.	Instep Circum.	Toe Depth
3	4	128	60	156-168	165-177	23-29
4	5	136	62	162-174	171-183	24-30
5	6	145	65	168-180	177-189	24-30
6	7	153	67	175-187	184-196	25-31
7	8	162	69	181-193	190-202	25-31
8	9	170	71	187-199	196-208	25-31
9	10	179	73	193-206	203-215	26-32
10	11	187	75	200-212	209-221	26-32
11	12	196	77	206-218	215-227	26-32
12	13	204	79	213-225	222-234	27-33
13	1	213	81	219-231	228-240	27-33
1	2	221	83	226-238	235-247	28-34
2	3	230	86	232-244	241-253	28-34
3	4	238	88	238-250	247-259	28-34
4	5	247	90	245-257	254-266	29-35
5	6	255	92	251-263	260-272	29-35
6	7	263	94	257-269	266-278	29-35
7	8	272	96	264-276	273-285	30-36
8	9	280	98	270-282	281-286	30-36
9	10	288	100	276-288	289-294	31-37

<b>Sizes available (UK)</b>	<b>Child</b>	<b>Adolescent</b>
<b>Grey Trainer</b>	<b>3-13</b>	<b>1-7</b>
<b>BlackTrainer</b>	<b>10-13</b>	<b>1-8</b>
<b>School Shoe</b>	<b>9-13</b>	<b>1-9</b>

**Please advise Child or Adolescent when ordering.**

# Bespoke Footwear

All Steprite styles can be made to bespoke specifications.

For addition of insoles please note that no allowance is built in for Insoles/TCI's/FFO's.

In the event that insoles/TCI's/FFO's have to be incorporated into the footwear and are NOT being made by our technicians, the allowance required for accommodation MUST be supplied with the footwear.

If outline of feet are of different lengths, please state whether the footwear is to be made the same length for both feet or if different lengths are required.

We are always happy to make footwear to casts if patients have difficult feet to draft. Please, however, supply measures and boot lengths if applicable.

Please order a plain upper style for the patients first pair so the correct fitting can be established.

I.e. Austin and Houston are problematic due to their complicated patterns.

Some additional styles more suitable for bespoke manufacture are shown below.



**Clyde**  
Three eyelet derby shoe with plain vamps.



**Kelvin**  
Three eyelet derby shoe with semi-brogue design.



**Tay**  
Five eyelet derby shoe with plain vamps, micro platform sole.



**Moriston**  
Derby boot with padded collars (depth of collar optional, 25mm -50mm).



**Tummel**  
Desert boot with optional turn down collar or padded collar.



**Devon**  
Two eyelet derby shoe with low entry and mock apron stitch design, two tone material option.



### Spey

Strap and buckle fastening, sandal style shoe, two tone material optional. 072 lightweight micro sole and heel.



### Esk

Golosh lorne boot with padded collar, no side seams.



### Ruthven

Open to toe derby boot with padded collar, depth variable 25mm - 50mm, velcro fastening.



### Rannoch

Derby boot with deep padded collar to cushion ankle bones, 072 lightweight microsole and heel.



### Yew

Low opening derby shoe with padded collar, sole as OTIS but with cutaway waist.



### Larch

Derby ankle boot with padded collar and twin velcro strap fasten, lightweight TCR commando sole unit.



### Ash

Low open derby shoe with padded collar, twin velcro strap fastening, platform micro sole unit cutaway at waist and heel with ribbed stick on sole on top.



### Suilven

Derby shoe with velcro strap fastening, padded collar, new lightweight micro deep ribbed sole, affording more flexibility and bounce.

# Sole Options for Modular Footwear Only

Our standard soling is a lightweight composite material, which is harder wearing than EVA, gives better grip and is flexible. Sole units can be interchangeable with those illustrated on individual footwear styles. If changed however, these are non returnable.

## Standard Sole Units



Thru Wedge



Owen

## Optional Sole Units



Hard Wearing Sole



Eva



Thru Rocker

## Linings



### ecoTex™ BAMBUS

**Bambus:** Made from Bamboo. It has natural anti bacterial properties using no chemicals, is antistatic, strong and durable and has a natural deodorising property.



**Go Guard:** Developed for the diabetic foot, with silver ions woven into the fabric to kill existing bacteria and prevent regeneration and reduce odour without the use of chemicals.

# Insole System

Buchanan Orthotics provides a full range of custom functional foot orthoses and total contact insoles with an extensive range of additions allowing insoles tailored to each individual.

Our highly skilled technicians and professional expertise provides a great service and product range. Our customer service team is also on hand to provide any support or information required.

CAD CAM technology is used to allow accurate and repeatable insoles manufactured according to the prescription given by the professional.

All insole styles are fully custom made to the cast or impression box sent in with a completed insole specification form. Any required modification can be undertaken by our technicians and then transformed to the milling machine to carve out the shell style required. Any additions and top covers will then be added to complete the desired Orthotics insole.

We understand that it is important to take into consideration each individual's needs, therefore there is an extensive range of materials available to ensure exact specification is received.



## TCI's, Plastic Orthoses & FFO's



### Total Contact Insoles

We manufacture total contact insoles to your own specification and in a choice of materials with varying degrees of density.

However, if no preference is given, our standard is:

- Base Layer: 3mm MED EVA
- Mid Layer: 3mm PINK PORON
- Top Cover: Leather

When ordering TCI's and footwear a negative cast or impression box should accompany the order to enable us to achieve a perfect fit.

### FFOs

We manufacture functional foot orthotics to your specifications from casts or impression boxes.

Our FFO's are available in a variety of materials, for example, EVA, carbon fibre or polypropylene.

Additions and modifications to our FFO's are made exactly to your specification.



### Plastic Orthoses

We manufacture the full range of plastic orthoses.

AFO, HINGED AFO, GROUND REACTION AFO, KAFO, HKAFO, SPINAL ORTHOSES.

Illustration shows paediatric AFO with transfers and window.

# Fabric Supports



**Handzon Back Support**  
 Easy application for patients with limited hand mobility. Soft curved front ensures comfortable wear and elasticated back gives a snug fit.  
 Also available in extra depth. Back can be boned and have additional steels or be pocketed to include a cloth covered hot/cold pad.



**Standard Lumbar Support (right)**  
 From stock or made to measure.

**Wrist Support (below left & centre)**  
 Beige or Black, Standard or Long.

**Child's Gaiter (below far right)**  
 Available in any size.



# Colours

Softie Leather	Tan	Grey	White	Navy	Brown	Burgundy	Camel	Black
	Bone	Light Beige	Mushroom	Taupe	Dark Brown	Sage	Tosca	Light Blue
	Soft Grain Brown	Soft Grain Tan	Soft Grain Black	Pink	Jeans			
Calf	Taupe	Burgundy	Dark Brown	Brown	Black	Tan	Grey	Navy
Suede	Black	Navy	Grey	Camel	Taupe	Light Beige	Burgundy	Sage
	Tan	Brown	Dark Brown					
Nubuck	Navy	Black	Burgundy	Grey	Taupe	Brown		
Heavy Grain	Brown	Black						
Stretch Leather				Brown	Black	Navy	Taupe	Leopard Skin
Neoprene	Burgundy	Light Beige	Brown	Navy	Black			
Metallic	Cappuccino	Crushed Patent Black	Antique Silver	Peach				
Raindrop		Black	Brown	Tan				
Tweed	Brown	Black						

Colours shown are not exact.

# Ladies Regular Lasts

Features: round toe, 2 x 3mm poron insoles & footbed

Standard

Size	Heel Width	Joint Width	Joint Circum. with Poron insoles & Footbed	Joint Circum. without Poron insoles & Footbed	Instep Circum. with Poron insoles & Footbed	Instep Circum. without Poron insoles & Footbed	Toe Depth	Last Length
1	57	85	217	235	201	219	24-33	223
2	58	88	224	242	208	226	25-34	231
3	60	90	230	248	214	232	26-35	239
4	61	93	237	255	221	239	26-35	248
5	63	95	243	261	227	245	27-36	257
6	64	98	250	268	234	252	28-37	265
7	66	100	256	274	240	258	29-38	274
8	67	103	263	281	247	265	29-38	282
9	69	105	269	287	253	271	30-39	291
10	70	108	276	294	260	278	31-40	300

Wide

Size	Heel Width	Joint Width	Joint Circum. with Poron insoles & Footbed	Joint Circum. without Poron insoles & Footbed	Instep Circum. with Poron insoles & Footbed	Instep Circum. without Poron insoles & Footbed	Toe Depth	Last Length
1	61	91	229	247	209	227	24-33	223
2	62	94	236	254	216	234	25-34	231
3	64	96	242	260	222	240	26-35	239
4	65	99	249	267	229	247	26-35	248
5	67	101	255	273	235	253	27-36	257
6	68	104	262	280	242	260	28-37	265
7	70	106	268	286	248	266	29-38	274
8	71	109	275	293	255	273	29-38	282
9	73	111	281	299	261	279	30-39	291
10	74	114	288	306	268	286	31-40	300

Slim

Size	Heel Width	Joint Width	Joint Circum. with Poron insoles & Footbed	Joint Circum. without Poron insoles & Footbed	Instep Circum. with Poron insoles & Footbed	Instep Circum. without Poron insoles & Footbed	Toe Depth	Last Length
1	53	78	205	223	193	212	24-33	223
2	54	82	212	230	200	218	25-34	231
3	56	84	218	236	206	224	26-35	239
4	57	87	225	243	213	231	26-35	248
5	59	89	231	249	219	237	27-36	257
6	60	92	238	256	226	244	28-37	265
7	62	94	244	262	232	250	29-38	274
8	63	97	251	269	239	257	29-38	282
9	65	99	257	275	245	263	30-39	291
10	66	102	264	282	252	270	31-40	300

Deep

Size	Heel Width	Joint Width	Joint Circum. with Poron insoles & Footbed	Joint Circum. without Poron insoles & Footbed	Instep Circum. with Poron insoles & Footbed	Instep Circum. without Poron insoles & Footbed	Toe Depth	Last Length
1	57	85	229	247	201	209	30-39	223
2	58	88	236	254	208	226	31-40	231
3	60	90	242	260	214	232	32-41	239
4	61	93	249	267	221	239	32-41	248
5	63	95	255	270	227	245	33-42	257
6	64	98	262	280	234	252	34-43	265
7	66	100	268	286	240	258	35-44	274
8	67	103	275	293	247	265	35-44	282
9	69	105	281	299	253	271	36-45	291
10	70	108	288	306	260	278	37-46	300

Wide & Deep

Size	Heel Width	Joint Width	Joint Circum. with Poron insoles & Footbed	Joint Circum. without Poron insoles & Footbed	Instep Circum. with Poron insoles & Footbed	Instep Circum. without Poron insoles & Footbed	Toe Depth	Last Length
1	61	91	241	259	209	227	30-39	223
2	62	94	248	266	216	234	31-40	231
3	64	96	254	272	222	240	32-41	239
4	65	99	261	279	229	247	32-41	248
5	67	101	267	285	235	253	33-42	257
6	68	104	274	292	242	260	34-43	265
7	70	106	280	298	248	266	35-44	274
8	71	109	287	305	255	273	35-44	282
9	73	111	293	311	261	279	36-45	291
10	74	114	300	318	268	286	37-46	300



## Donna (Narrower Heel)

Size	Heel Width	Instep Circum. with Poron insoles & Footbed	Instep Circum. without Poron insoles & Footbed
1	-	-	-
2	54	198	217
3	56	204	223
4	57	211	230
5	59	217	236
6	60	224	243
7	62	230	249
8	63	237	256
9	-	-	-
10	-	-	-

Size	Heel Width	Instep Circum. with Poron insoles & Footbed	Instep Circum. without Poron insoles & Footbed
1	-	-	-
2	58	206	225
3	60	212	231
4	61	219	238
5	63	225	244
6	64	232	251
7	66	238	257
8	67	245	264
9	-	-	-
10	-	-	-

Size	Heel Width	Instep Circum. with Poron insoles & Footbed	Instep Circum. without Poron insoles & Footbed
1	-	-	-
2	50	190	209
3	52	196	215
4	53	203	222
5	55	209	228
6	56	216	235
7	58	222	241
8	59	229	248
9	-	-	-
10	-	-	-

Standard

Wide

Slim

## Diana (Shallower Forepart)

Size	Joint Circum. with Poron insoles & Footbed	Joint Circum. without Poron insoles & Footbed	Toe Depth
1	-	-	-
2	214	233	21-30
3	220	239	22-31
4	227	246	22-31
5	233	252	23-32
6	240	259	24-33
7	246	265	25-34
8	253	272	25-34
9	-	-	-
10	-	-	-

Standard

## Denise (Narrower Heel & Shallower Forepart)

Size	Heel Width	Joint Circum. with Poron insoles & Footbed	Joint Circum. without Poron insoles & Footbed	Instep Circum. with Poron insoles & Footbed	Instep Circum. without Poron insoles & Footbed	Toe Depth
1	-	-	-	-	-	-
2	54	214	233	198	217	21-30
3	56	220	239	204	223	22-31
4	57	227	246	211	230	22-31
5	59	233	252	217	236	23-32
6	60	240	259	224	243	24-33
7	62	246	265	230	249	25-34
8	63	253	272	237	256	25-34
9	-	-	-	-	-	-
10	-	-	-	-	-	-

Size	Heel Width	Joint Circum. with Poron insoles & Footbed	Joint Circum. without Poron insoles & Footbed	Instep Circum. with Poron insoles & Footbed	Instep Circum. without Poron insoles & Footbed	Toe Depth
1	-	-	-	-	-	27-36
2	58	226	245	206	225	27-36
3	60	232	251	212	231	28-37
4	61	239	258	219	238	29-38
5	63	245	264	225	244	29-38
6	64	252	271	232	251	30-39
7	66	258	277	238	257	31-40
8	67	265	284	245	264	31-40
9	-	-	-	-	-	32-41
10	-	-	-	-	-	-

Size	Joint Circum. without Poron insoles & Footbed	Joint Circum. with Poron insoles & Footbed	Toe Depth
1	-	-	-
2	226	245	21-30
3	232	251	22-31
4	239	258	22-31
5	245	264	23-32
6	252	271	24-33
7	258	277	25-34
8	265	284	25-34
9	-	-	-
10	-	-	-

Wide

Size	Joint Circum. with Poron insoles & Footbed	Joint Circum. without Poron insoles & Footbed	Toe Depth
1	-	-	-
2	202	221	21-30
3	208	227	22-31
4	215	234	22-31
5	221	240	23-32
6	228	247	24-33
7	234	253	25-34
8	241	260	25-34
9	-	-	-
10	-	-	-

Slim

Size	Heel Width	Joint Circum. with Poron insoles & Footbed	Joint Circum. without Poron insoles & Footbed	Instep Circum. with Poron insoles & Footbed	Instep Circum. without Poron insoles & Footbed	Toe Depth
1	-	-	-	-	-	-
2	50	202	221	190	209	21-30
3	52	208	227	196	215	22-31
4	53	215	234	203	222	22-31
5	55	221	240	209	228	23-32
6	56	228	247	216	235	24-33
7	58	234	253	222	241	25-34
8	59	241	260	229	248	25-34
9	-	-	-	-	-	-
10	-	-	-	-	-	-

# Regular Lasts

Features: round toe, 2 x 3mm poron insoles & footbed

Standard

Size	Heel Width	Joint Width	Joint Circum. with Poron Insoles & Footbed	Joint Circum. without Poron insoles & Footbed	Instep Circum. with Poron insoles & Footbed	Instep Circum. without Poron insoles & Footbed	Toe Depth	Last Length
4	62	91	231	249	221	239	31-40	246
5	64	93	237	255	227	245	31-40	254
6	65	96	244	262	234	252	32-41	263
7	67	98	250	268	240	258	33-42	271
8	68	101	257	275	247	265	33-42	279
9	70	103	263	281	253	271	34-43	288
10	71	106	270	288	260	278	35-44	296
11	73	108	276	294	266	284	35-44	305
12	74	111	283	301	273	291	36-45	314
13	76	113	289	307	279	297	37-46	322
14	77	116	296	314	286	304	37-46	331

Wide

Size	Heel Width	Joint Width	Joint Circum. with Poron Insoles & Footbed	Joint Circum. without Poron insoles & Footbed	Instep Circum. with Poron insoles & Footbed	Instep Circum. without Poron insoles & Footbed	Toe Depth	Last Length
4	66	97	243	261	229	247	31-40	246
5	68	99	249	267	235	253	31-40	254
6	69	102	256	274	242	260	32-41	263
7	71	104	262	280	248	266	33-42	271
8	72	107	269	287	255	273	33-42	279
9	74	109	275	293	261	279	34-43	288
10	75	112	282	300	268	286	35-44	296
11	77	114	288	306	274	292	35-44	305
12	78	117	295	313	281	299	36-45	314
13	80	119	301	319	287	305	37-46	322
14	81	122	308	326	294	312	37-46	331

Slim

Size	Heel Width	Joint Width	Joint Circum. with Poron Insoles & Footbed	Joint Circum. without Poron insoles & Footbed	Instep Circum. with Poron insoles & Footbed	Instep Circum. without Poron insoles & Footbed	Toe Depth	Last Length
4	58	85	219	237	213	231	31-40	246
5	60	87	225	243	219	237	31-40	254
6	61	90	232	250	226	244	32-41	263
7	63	92	238	256	232	250	33-42	271
8	64	95	245	263	239	257	33-42	279
9	66	97	251	269	245	263	34-43	288
10	67	100	258	276	252	270	35-44	296
11	69	102	264	282	258	276	35-44	305
12	70	105	271	289	265	283	36-45	314
13	72	107	277	295	271	289	37-46	322
14	73	110	284	302	278	296	37-46	331

Deep

Size	Heel Width	Joint Width	Joint Circum. with Poron Insoles & Footbed	Joint Circum. without Poron insoles & Footbed	Instep Circum. with Poron insoles & Footbed	Instep Circum. without Poron insoles & Footbed	Toe Depth	Last Length
4	62	91	243	261	221	239	37-46	246
5	64	93	249	267	227	245	37-46	254
6	65	96	256	274	234	252	38-47	263
7	67	98	262	280	240	258	39-48	271
8	68	101	269	287	247	265	39-48	279
9	70	103	275	293	253	271	40-49	288
10	71	106	282	300	260	278	41-50	296
11	73	108	288	306	266	284	41-50	305
12	74	111	295	313	273	291	42-51	314
13	76	113	301	319	279	297	43-52	322
14	77	116	308	326	286	304	43-52	331

Wide & Deep

Size	Heel Width	Joint Width	Joint Circum. with Poron Insoles & Footbed	Joint Circum. without Poron insoles & Footbed	Instep Circum. with Poron insoles & Footbed	Instep Circum. without Poron insoles & Footbed	Toe Depth	Last Length
4	66	97	255	273	229	247	37-46	246
5	68	99	261	279	235	253	37-46	254
6	69	102	268	286	242	260	38-47	263
7	71	104	274	292	248	266	39-48	271
8	72	107	281	299	255	273	39-48	279
9	74	109	287	305	261	279	40-49	288
10	75	112	294	312	268	286	41-50	296
11	77	114	300	318	274	292	41-50	305
12	78	117	307	325	281	299	42-51	314
13	80	119	313	331	287	305	43-52	322
14	81	122	320	338	294	312	43-52	331

### Albert (Narrower Heel)

Size	Heel Width	Instep Circum. with Poron insoles & Footbed	Instep Circum. without Poron insoles & Footbed
4	58	211	230
5	60	217	236
6	61	224	243
7	63	230	249
8	64	237	256
9	66	243	262
10	67	250	269
11	69	256	275
12	70	263	282
13	-	-	-
14	-	-	-

### Wide

Size	Heel Width	Instep Circum. with Poron insoles & Footbed	Instep Circum. without Poron insoles & Footbed
4	62	219	238
5	64	225	244
6	65	232	251
7	67	238	257
8	68	245	264
9	70	251	270
10	71	258	277
11	73	264	283
12	74	271	290
13	-	-	-
14	-	-	-

### Slim

Size	Heel Width	Instep Circum. with Poron insoles & Footbed	Instep Circum. without Poron insoles & Footbed
4	58	211	230
5	60	217	236
6	61	224	243
7	63	230	249
8	64	237	256
9	66	243	262
10	67	250	269
11	69	256	275
12	70	263	282
13	-	-	-
14	-	-	-

### George (Shallower Forepart)

Size	Joint Circum. with Poron insoles & Footbed	Joint Circum. without Poron insoles & Footbed	Toe Depth
4	221	240	27-36
5	227	246	27-36
6	234	253	28-37
7	240	259	29-38
8	247	266	29-38
9	253	272	30-39
10	260	279	31-40
11	266	285	31-40
12	273	292	32-41
13	-	-	-
14	-	-	-

### Richard (Narrower Heel & Shallower Forepart)

Size	Heel Width	Joint Circum. with Poron insoles & Footbed	Joint Circum. without Poron insoles & Footbed	Instep Circum. with Poron insoles & Footbed	Instep Circum. without Poron insoles & Footbed	Toe Depth
4	58	221	240	211	230	27-36
5	60	227	246	217	236	27-36
6	61	234	253	224	243	28-37
7	63	240	259	230	249	29-38
8	64	247	266	237	256	29-38
9	66	253	272	243	262	30-39
10	67	260	279	250	269	31-40
11	69	266	285	256	275	31-40
12	70	273	292	263	282	32-41
13	-	-	-	-	-	-
14	-	-	-	-	-	-

### Wide

Size	Joint Circum. with Poron insoles & Footbed	Joint Circum. without Poron insoles & Footbed	Toe Depth
4	233	252	27-36
5	239	258	27-36
6	246	265	28-37
7	252	271	29-38
8	259	278	29-38
9	265	284	30-39
10	272	291	31-40
11	278	297	31-40
12	285	304	32-41
13	-	-	-
14	-	-	-

### Slim

Size	Heel Width	Joint Circum. with Poron insoles & Footbed	Joint Circum. without Poron insoles & Footbed	Instep Circum. with Poron insoles & Footbed	Instep Circum. without Poron insoles & Footbed	Toe Depth
4	62	233	252	219	236	27-36
5	64	239	258	225	244	27-36
6	65	246	265	232	251	28-37
7	67	252	271	238	257	29-38
8	68	259	278	245	264	29-38
9	70	265	284	251	270	30-39
10	71	272	291	258	277	31-40
11	73	278	297	264	283	31-40
12	74	285	304	271	290	32-41
13	-	-	-	-	-	-
14	-	-	-	-	-	-

### Slim

Size	Joint Circum. with Poron insoles & Footbed	Joint Circum. without Poron insoles & Footbed	Toe Depth
4	209	228	27-36
5	215	234	27-36
6	222	241	28-37
7	228	247	29-38
8	235	254	29-38
9	241	260	30-39
10	248	267	31-40
11	254	273	31-40
12	261	280	32-41
13	-	-	-
14	-	-	-

### Slim

Size	Heel Width	Joint Circum. with Poron insoles & Footbed	Joint Circum. without Poron insoles & Footbed	Instep Circum. with Poron insoles & Footbed	Instep Circum. without Poron insoles & Footbed	Toe Depth
4	54	209	228	203	222	27-36
5	56	215	234	209	228	27-36
6	57	222	241	216	235	28-37
7	59	228	247	222	241	29-38
8	60	235	254	229	248	29-38
9	62	241	260	235	254	30-39
10	63	248	267	242	261	31-40
11	65	254	273	248	267	31-40
12	66	261	280	255	274	32-41
13	-	-	-	-	-	-
14	-	-	-	-	-	-

# Ladies Contour Lasts

## Regular Contour Lasts

Features: round toe, 2 x 3mm poron insoles & footbed

Standard

Size	Heel Width	Joint Width	Joint Circum. with Poron insoles & Footbed	Joint Circum. without Poron insoles & Footbed	Instep Circum. with Poron insoles & Footbed	Instep Circum. without Poron insoles & Footbed	Toe Depth	Last Length
1	55	85	212	224	206	218	24-30	222
2	57	87	218	230	212	224	25-31	230
3	59	89	224	235	218	230	26-32	238
4	61	92	230	242	224	236	26-32	247
5	63	94	237	249	231	243	27-33	256
6	65	95	243	255	237	249	28-34	265
7	67	97	249	261	243	255	29-35	274
8	68	100	255	267	250	262	30-36	282
9	70	102	260	272	256	268	31-37	291
10	71	105	267	279	262	274	31-37	300
11	73	107	273	285	268	280	33-39	309
12	75	110	280	292	275	287	34-40	318
13	76	112	291	303	283	295	35-41	326

## Donna Contour (Narrower Heel)

Size	Heel Width	Instep Circum. with Poron insoles & Footbed	Instep Circum. without Poron insoles & Footbed
1	51	203	215
2	53	210	222
3	55	216	228
4	57	221	233
5	59	227	239
6	61	232	244
7	63	237	249
8	64	242	254
9	66	248	260
10	67	255	267
11	69	262	274
12	71	268	280
13	72	273	285

Wide

Size	Heel Width	Joint Width	Joint Circum. with Poron Insoles & Footbed	Joint Circum. without Poron insoles & Footbed	Instep Circum. with Poron insoles & Footbed	Instep Circum. without Poron insoles & Footbed	Toe Depth	Last Length
1	61	91	229	241	219	231	24-30	222
2	63	94	235	247	225	237	25-31	230
3	65	96	241	253	231	243	26-32	238
4	66	99	247	259	237	249	26-32	247
5	68	101	253	265	243	255	27-33	256
6	69	103	259	271	249	261	28-34	265
7	71	105	265	277	255	267	29-35	274
8	72	107	271	283	261	273	30-36	282
9	74	110	277	289	267	279	31-37	291
10	75	113	284	296	274	286	31-37	300
11	77	115	290	302	280	292	33-39	309
12	78	118	296	308	287	299	34-40	318
13	80	120	302	314	293	305	35-41	326

Size	Heel Width	Instep Circum. with Poron insoles & Footbed	Instep Circum. without Poron insoles & Footbed
1	57	215	227
2	59	219	231
3	61	225	237
4	62	230	242
5	64	236	248
6	65	242	254
7	67	248	260
8	68	254	266
9	70	260	272
10	71	266	278
11	73	273	285
12	74	279	291
13	76	286	298

Slim

Size	Heel Width	Joint Width	Joint Circum. with Poron Insoles & Footbed	Joint Circum. without Poron Insoles & Footbed	Instep Circum. with Poron insoles & Footbed	Instep Circum. without Poron insoles & Footbed	Toe Depth	Last Length
1	53	79	204	216	199	211	24-30	222
2	54	82	210	222	205	217	25-31	230
3	56	85	216	228	211	223	26-32	238
4	57	87	222	234	217	229	26-32	247
5	59	89	228	240	223	235	27-33	256
6	60	92	234	246	229	241	28-34	265
7	62	94	240	252	235	247	29-35	274
8	63	97	246	258	241	253	30-36	282
9	65	99	252	264	247	259	31-37	291
10	66	102	258	270	253	265	31-37	300
11	68	105	264	276	259	271	33-39	309
12	70	107	270	282	265	276	34-40	318
13	71	109	276	288	271	283	35-41	326

Size	Heel Width	Instep Circum. with Poron insoles & Footbed	Instep Circum. without Poron insoles & Footbed
1	49	193	205
2	50	199	211
3	52	206	218
4	53	212	224
5	55	218	230
6	56	224	236
7	58	230	242
8	59	236	248
9	61	242	254
10	62	248	260
11	64	254	266
12	66	260	272
13	67	266	279

# Gents Contour Lasts

Standard

## Regular Contour Lasts

Features: round toe, 2 x 3mm poron insoles & footbed

Size	Heel Width	Joint Width	Joint Circum. with Poron insoles & Footbed	Joint Circum. without Poron insoles & Footbed	Instep Circum. with Poron insoles & Footbed	Instep Circum. without Poron insoles & Footbed	Toe Depth	Last Length
1	55	85	212	225	206	218	24-30	222
2	57	87	218	230	212	224	25-31	230
3	59	89	224	235	218	230	26-32	238
4	61	92	230	242	224	236	26-32	247
5	63	94	237	249	231	243	27-33	256
6	65	95	243	255	237	249	28-34	265
7	67	97	249	261	243	255	29-35	274
8	68	100	255	267	250	262	30-36	282
9	70	102	260	272	256	268	31-37	291
10	71	105	267	279	262	274	31-37	300
11	73	107	273	285	268	280	33-39	309
12	75	110	280	292	275	287	34-40	318
13	76	112	291	303	283	295	35-41	326
14	-	-	-	-	-	-	-	-

## Albert Contour (Narrower Heel)

Size	Heel Width	Instep Circum. with Poron insoles & Footbed	Instep Circum. without Poron insoles & Footbed
1	51	203	215
2	53	210	222
3	55	216	228
4	57	221	233
5	59	227	239
6	61	232	244
7	63	237	249
8	64	242	254
9	66	248	260
10	67	255	267
11	69	262	274
12	71	268	280
13	72	273	285
14	-	-	-

Wide

Size	Heel Width	Joint Width	Joint Circum. with Poron Insoles & Footbed	Joint Circum. without Poron Insoles & Footbed	Instep Circum. with Poron insoles & Footbed	Instep Circum. without Poron insoles & Footbed	Toe Depth	Last Length
1	61	91	229	241	213	219	24-30	222
2	63	94	235	247	225	237	25-31	230
3	65	96	241	253	231	243	26-32	238
4	66	99	247	259	237	249	26-32	247
5	68	101	253	265	243	255	27-33	256
6	69	103	259	271	249	261	28-34	265
7	71	105	265	277	255	267	29-35	274
8	72	107	271	283	261	273	30-36	282
9	74	110	277	289	267	279	31-37	291
10	75	113	284	296	274	286	31-37	300
11	77	115	290	302	280	292	33-39	309
12	78	118	296	308	287	299	34-40	318
13	80	120	302	314	293	305	35-41	326
14	-	-	-	-	-	-	-	-

Size	Heel Width	Instep Circum. with Poron insoles & Footbed	Instep Circum. without Poron insoles & Footbed
1	57	215	227
2	59	219	231
3	61	225	237
4	62	230	242
5	64	236	248
6	65	242	254
7	67	248	260
8	68	254	266
9	70	260	272
10	71	266	278
11	73	273	285
12	74	279	291
13	76	286	298
14	-	-	-

Slim

Size	Heel Width	Joint Width	Joint Circum. with Poron insoles & Footbed	Joint Circum. without Poron insoles & Footbed	Instep Circum. with Poron insoles & Footbed	Instep Circum. without Poron insoles & Footbed	Toe Depth	Last Length
1	53	79	204	216	199	211	24-30	222
2	54	82	210	222	205	217	25-31	230
3	56	85	216	228	211	223	26-32	238
4	57	87	222	234	217	229	26-32	247
5	59	89	228	240	223	235	27-33	256
6	60	92	234	246	229	241	28-34	265
7	62	94	240	252	235	247	29-35	274
8	63	97	246	258	241	253	30-36	282
9	65	99	252	264	247	259	31-37	291
10	66	102	258	270	253	265	31-37	300
11	68	105	264	276	259	271	33-39	309
12	70	107	270	282	265	276	34-40	318
13	71	109	276	288	271	283	35-41	326
14	-	-	-	-	-	-	-	-

Size	Heel Width	Instep Circum. with Poron insoles & Footbed	Instep Circum. without Poron insoles & Footbed
1	49	193	205
2	50	199	211
3	52	206	218
4	53	212	224
5	55	218	230
6	56	224	236
7	58	230	242
8	59	236	248
9	61	242	254
10	62	248	260
11	64	254	266
12	66	260	272
13	67	266	279
14	-	-	-

# Last Adaptions

**Half  
Sizes Also  
Available**

**Maximum Additional  
Instep Circumference**  
All sizes 50mm extra  
to catalogue.

**Joint Width**  
Where additional Joint Width  
above 4mm is required and  
unless otherwise advised,  
we will automatically  
manufacture footwear with  
bespoke soles and this will  
incur an extra charge.

**Maximum Additional  
Forepart Depth**  
Regular last: 10mm.  
Contour last: 6mm.

**Bunion Pocket**  
Maximum 9mm in  
3mm increments.

**Standard Toe Depth**  
Sizes 1-3: additional 10mm.  
Size 4 and upward: maximum 52mm in total.

**Contour Toe Depth**  
All sizes: additional 6mm.

**Taylor Bunion (bunionette) Pocket**  
Maximum 6mm in 3mm increments.

**Quarter Height**  
Can be lowered at the ankle.  
Give height from floor to base of  
ankle or specify quarter height.

**Standard Toe Shape**

**Contour Toe Shape**

**Standard Toe Spring**

**Contour Toe Spring**

**Rough Fit**  
Additional individual modifications after rough fit can be undertaken up to a maximum of 6mm, ie. instep, toe box, length etc., though not a combination of these.

**Major Upper Modifications**  
Must be requested at the ordering stage.

Gents UK Sizes										
5	5½	6½	7	8	9	9½	10½	11	11½	12
Gents European Sizes										
38	39	40	41	42	43	44	45	46	47	48

Ladies UK Sizes									
3	4	5	5½	6½	7	8	9	9½	
Ladies European Sizes									
36	37	38	39	40	41	42	43	44	

## Modular Footwear Conditions of Sale

Available at no extra cost

- Half Sizes
- Touch & Close Fastening
- Trial Fitting Service

### Despatch

Modular footwear will be despatched within seven to ten working days. Shoes requiring adaptations will be despatched within ten to fifteen working days.

As an orthotic company we know the importance of meeting our delivery dates and are always striving to improve on them. On occasion, however, circumstances beyond our control may result in a slight delay and our customer support team will contact you with a revised despatch date should a delay occur.

### Carriage

Carriage free for Northern Ireland and UK mainland.

### Payment

Payment is due thirty days from invoice date if account facilities are established.

Volume discounts are available for orders of 20, 40 and 80 pairs per month.

20 pairs per month - 5% discount  
40 pairs per month - 10% discount  
80 pairs per month - 15% discount

### Returns

All shoe styles on regular lasts without adaptations/modifications are returnable for credit or exchange if returned within two months in original packaging and in full resaleable condition.

Boot styles are subject to a 30% handling charge.

Split sizes or fittings are non returnable.

Any footwear requested with a through sole will not be returnable unless the through sole is standard for that type of footwear.

### Ownership of Goods

The footwear shall remain the absolute property of Buchanan Orthotics Limited, until such time as the customer shall have paid the agreed price in full.

### Colours

Our colour swatches shown here may vary from the actual shades due to the limitations of the reproduction processes.

Brochure design and art direction by Source Design, Glasgow. Tel 0141 353 5320 Photography by Blackdoor Studio Tel 07503 664 823



**Buchanan Orthotics**

Buchanan Orthotics Ltd.  
603 Helen Street, Glasgow, G51 3AR.  
Tel: 0141 440 1999.  
Fax: 0141 440 5999.  
Email: [sales@buchananorthotics.co.uk](mailto:sales@buchananorthotics.co.uk)  
Web: [www.buchananorthotics.co.uk](http://www.buchananorthotics.co.uk)



Certificate no. FM 55047

